



Recombinant Human Proline-rich protein 23C (PRR23C)

Product Code	CSB-BP754559HU
Abbreviation	PRR23C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZRP0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGSRPCSPSA CLAPWWGQQP GGP GPAKR SR LEEPAGPESR AAPSPEDPAG TPAVDALTSM VVLDAGCALR VPLEDVDLVL ELAPMSVLRV SLGGHTLIVI PEVLLSSVDE CSGAQGDWSA GLEVDVFLGA HGEDVVVEQE VCASVPEIAA EEEAYEEDAD SEFPELWMDS AAGSAAGLYP SARSMFSPYR EGPIRGPCAL APNPSSERRS PRPIFDLEFH LLEPVPSSPL QPLPPSPSPG PHARPELPER PPCKVRRRLF QE
Source	Baculovirus
Target Names	PRR23C
Protein Names	Recommended name: Proline-rich protein 23C
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.