



Recombinant Human Cyclic AMP-responsive element-binding protein 3-like protein 2 (CREB3L2), partial, Yeast

Product Code	CSB-YP754617HU
Abbreviation	CREB3L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q70SY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	CREB3L2
Protein Names	Recommended name: Cyclic AMP-responsive element-binding protein 3-like protein 2 Short name= cAMP-responsive element-binding protein 3-like protein 2 Alternative name(s): BBF2 human homolog on chromosome 7 Cleaved into the following chain: 1.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	full length of BC063666
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.