



Recombinant Human immunodeficiency virus type 2 subtype B Protein Rev (rev)

Product Code	CSB-MP754993HME
Abbreviation	rev
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q74232
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 2 subtype B (isolate EHO) (HIV-2)
Purity	≥85% (SDS-PAGE)
Sequence	MNARERDLQK GLRLLHLLHQ TNPYPQGP GT ASQRRNR RRR WKQRGLQILA LADRIHPLPD SPTEGPLDLA IQRLQNLIK DLPNPPTSTP TAQASTCIPP IWDQLVPRSN PSSSQGCGRD SCERGEDLVG SPQESGRRDH CNTQEDQTRG
Source	Mammalian cell
Target Names	rev
Protein Names	Recommended name: Protein Rev Alternative name(s): Regulator of expression of viral proteins
Expression Region	1-150
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.