



Recombinant Human E3 ubiquitin-protein ligase Hakai (CBLL1)

Product Code	CSB-YP755210HU
Abbreviation	CBLL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75N03
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDHTDNELQG TNSSGSLGGL DVRRRIPIKL ISKQANKAKP APRTQRTINR MPAKAPPGDE EGFYNEEER YDCKGGELFA NQRRFPGHLF WDFQINILGE KDDTPVHFCD KCGLPKIYGM RMIPCKHVFC YDCAILHEKK GDKMCPGCSD PVQRIEQCTR GSLFMCSIVQ GCKRXYLSQR DLQAHINHRH MRAGKPVTRA SLENVHPPIA PPPTEIPERF IMPPDKHHMS HIPPKQHIMM PPPPLQHVPH EHYNQPHEDI RAPPAAELSM PPPPRSVSQE TFRISTRKHS NLITVPIQDD SNSGAREPPP PAPAPAHHP EYQGQPVVSH PHHIMPPQQH YAPPPPPPPP ISHPMPPHPQ AAGTPHLVYS QAPPPMPTSA PPPITPPPGH IIAQMPPYMN HPPPGPPPPQ HGGPPVTAPP PHHYNPNSLP QFTEDQGTLSP PFTQPGGMS PGIWPAPRGP PPPRLQGPP SQTPLPGPHH PDQTRYRPPY Q
Source	Yeast
Target Names	CBLL1
Protein Names	Recommended name: E3 ubiquitin-protein ligase Hakai EC= 6.3.2.- Alternative name(s): Casitas B-lineage lymphoma-transforming sequence-like protein 1 RING finger protein 188 c-Cbl-like protein 1
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.