



Recombinant Human immunodeficiency virus type 2 subtype B Virion infectivity factor (vif)

Product Code	CSB-BP755229HDAI
Abbreviation	vif
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q76635
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 2 subtype B (isolate UC1) (HIV-2)
Purity	>85% (SDS-PAGE)
Sequence	MEGEKNWIVV PTWRIPGRLE KWHSLVKYLK HRTKELQQVS YVPHHKVGVWA WWTCSRVIFFP LKEEAYLEVQ GYWNLTPERG FLSSYAVRLT WYKRSFYTDV TPDVADQLLH GSYFSCFTAN EVRRAIRGEK ILSYCNYP SA HEGQVPSLQF LALRVIQEGK DGSQGESATR KQRRRNRRS IRLARKNNNR AQQGSSQPLA PRTHFPGLAE VLGILA
Source	Baculovirus
Target Names	vif
Protein Names	Recommended name: Virion infectivity factor Short name= Vif Alternative name(s): Q protein SOR protein
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.