



Recombinant Human Protein FAM5C (FAM5C)

Product Code	CSB-EP755233HU-B
Abbreviation	FAM5C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q76B58
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AVSDQHA TSPFDWLLSD KGPFHRSQEY TDFVDRSRQG FSTRYKIYRE FGRWKVNNLA VERRNFLGSP LPLAPEFFRN IRLLGRRPTL QQITENLIKK YGTHFLLSAT LGGEESLTIF VDKRKLSKRA EGSDSTTNSS SVTLETLHQL AASYFIDRDS TLRRLLHHIQI ASTAIKV TET RTGPLGCSNY DNLDVSSVL VQSPENKIQL QGLQVLLPDY LQERFVQAAL SYIACNSEGE FICKENDCWC HCGPKFPECN CPSMDIQAME ENLLRITETW KAYNSDFEES DEFKLFMKRL PMNYFLNTST IMHLWTMDSN FQRRYELEN SMKQLFLKAQ KIVHKLFSL KRCHKQPLIS LPRQRTSTYW LTRIQSFLYC NENGLLGSFS EETHSCTCPN DQVVCTAFLP CTVG DASACL TCAPDNRTRC GTCNTGYMLS QGLCKPEVAE STDHYIGFET DLQDLEMKYL LQKTDRRIEV HAIFISNDRM LNSWFDPSWR KRMLLTLKSN KYKSSLVHMI LGLSLQICLT KNSTLEPVLA VYVNPFGGSH SESWFMPVNE NSFPDWERTK LDLPLQCYNW TLTLGNKWKT FFETVHIYLR SRIKSNPNG NESIYYEPL FIDPSRNLGY MKINNIQVFG YSMHFDPEAI RDLILQLDYP YTQGSQDSAL LQLLEIRDRV NKLSPPGQRR LDLFSCLLRH RLKLSTSEVV RIQSALQAFN AKLPNTMDYD TTKLCS
Source	E.coli
Target Names	BRINP3
Protein Names	Recommended name: Protein FAM5C Alternative name(s): DBCCR1-like protein 1
Expression Region	34-766
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.