



# Recombinant Human Protein FAM5C (FAM5C)

<b>Product Code</b>	CSB-MP755233HU
<b>Abbreviation</b>	FAM5C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q76B58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AVSDQHA TSPFDWLLSD KGPFHRSQEY TDFVDRSRQG FSTRYKIYRE FGRWKVNNLA VERRNFLGSP LPLAPEFFRN IRLLGRRPTL QQITENLIKK YGTHFLLSAT LGGEESLTIF VDKRKLSKRA EGSDSTTNSS SVTLETLHQL AASYFIDRDS TLRLHHIQI ASTAIKVTTET RTGPLGCSNY DNLDVSSVL VQSPENKIQL QGLQVLLPDY LQERFVQAAL SYIACNSEGE FICKENDCWC HCGPKFPECN CPSMDIQAME ENLLRITETW KAYNSDFEES DEFKLFMKRL PMNYFLNTST IMHLWTMDSN FQRRYELEN SMKQLFLKAQ KIVHKLFSL KRCHKQPLIS LPRQRTSTYW LTRIQSFLYC NENGLLGSFS EETHSCTCPN DQVVCTAFLP CTVGDAACL TCAPDNRTRC GTCNTGYMLS QGLCKPEVAE STDHYIGFET DLQDLEMKYL LQKTDRRIEV HAIFISNDRM LNSWFDPSWR KRMLLTLKSN KYKSSLVHMI LGLSLQICLT KNSTLEPVLA VYVNPFGGSH SESWFMPVNE NSFPDWERTK LDLPLQCYNW TLTLGNKWKT FFETVHIYLR SRIKSNPNG NESIYYEPL FIDPSRNLGY MKINNIQVFG YSMHFDPEAI RDLILQLDYP YTQGSQDSAL LQLLEIRDRV NKLSPPGQRR LDLFSCLLRH RLKLSTSEVV RIQSALQAFN AKLPNTMDYD TTKLCS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BRINP3
<b>Protein Names</b>	Recommended name: Protein FAM5C Alternative name(s): DBCCR1-like protein 1
<b>Expression Region</b>	34-766
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.