



Recombinant *Drosophila melanogaster* 39S ribosomal protein L41, mitochondrial (mRpL41)

Product Code	CSB-EP755468DLU
Abbreviation	mRpL41
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7JZM8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	GKRN FRKFNAYNKR GTRVVKEAQK TLANPPVAIH KRGVVRTGIL VDGQYVEIPE KIPDIIVPDL TGCKLKPYVS YKAPDVVQSE FTSLDLFNAV YSQKIIEDFK AGKLQKDGSA KEPSVNEQLT PEEALQRARK TGSDIF
Source	E.coli
Target Names	mRpL41
Protein Names	Recommended name: 39S ribosomal protein L41, mitochondrial Short name= L41mt Short name= MRP-L41
Expression Region	27-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.