



Recombinant Human Rab11 family-interacting protein 2 (RAB11FIP2)

Product Code	CSB-MP755487HU
Abbreviation	RAB11FIP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L804
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMLSEQAQKW FPTHVQVTVL QAKDLKPKGK SGTNDTYTII QLGKEKYSTS VAEKTLEPVW KEEASFELPG LLIQGSPEKY ILFLIVMHRS LVGLDKFLGQ VAINLNDIFE DKQRRKTEWF RLESKQGKRI KNRGEIKVNI QFMRNNMTAS MFDLSMKDKT RSPFAKLKDK MKGRKNDGTF SDTSSAIIPS THMPDANSEF SSGEIQMKSK PKKPFLGQP RLSSAHMSMSD LSGSHMSSEK LKAGTIGQTH LLGHQLDSFG TVPESGSLKS PHRRTLSFDT SKMNQPDSIV DEGELCFGRQ NDPFTNVTAS LPQKFATLPR KKNPFEESE TWDSSMNLFS KPIEIRKENK REKREKVSFLF ERVTGKKDSR RSDKLNNNGS DSPCDLKSPN AFSENQRQDYF DYESTNPFTA KFRASNIMPS SSFHMSPTSN EDLRKIPDSN PFDATAGYRS LTYEEVLQEL VKHKELLRRK DTHIRELEDY IDNLLVRVME ETPSILRVPY EPSRKAGKFS NS
Source	Mammalian cell
Target Names	RAB11FIP2
Protein Names	Recommended name: Rab11 family-interacting protein 2 Short name= Rab11-FIP2 Alternative name(s): NRip11
Expression Region	1-512
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.