



Recombinant Human SH2 domain-containing adapter protein F (SHF)

Product Code	CSB-MP755509HU
Abbreviation	SHF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7M4L6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQQEGGPVRS APCRTGTLEG SRQGSPGHRK RASPKGSLSS AQPHSWMLTP SPLNSHCAHR EPISSSPQPV ANGPKQKKKS NWRSTTRLRI IRLRDRLEPR PLAILEDYAD PFDVQETGEG SAGASGAPEK VPENDGYMEP YEAQKMMAEI RGSKETATQP LPLYDTPYEP EEDGATAEGE GAPWPRESRL PEDDERPPEE YDQPWEWKKE RISKAFAVDI KVIKDLPWPP PVGQLDSSPS LPDGDRDISG PASPLPEPSL EDSSAQFEGP EKSCLSPPGRE EKGRLLPRLS AGNPKSAKPL SMEPSSPLGE WTDPALPLEN QVWYHGAISR TDAENLLRLC KEASYLVRNS ETSKNDFSLK LKSSQGFMMH KLSRTKEHKY VLGQNSPPFS SVPEIVHHYA SRKLPIKGAE HMSLLYPVAI RTL
Source	Mammalian cell
Target Names	SHF
Protein Names	Recommended name: SH2 domain-containing adapter protein F
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.