



# Recombinant Human Putative ubiquitin carboxyl-terminal hydrolase 17-like protein 1 (USP17L1P)

<b>Product Code</b>	CSB-EP755915HU
<b>Abbreviation</b>	USP17L1P
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7RTZ2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGDDSLYLGG EWQFNHFSKL TSSRPDAAFA EIQRTSLPEK SPLSSETRVD LCDDLAPVAR QLAPREKLPL SSRPAAVGA GLQNMGNTCY ENASLQCLTY TLPLANYMLS REHSQTCQRP KCCMLCTMQA HITWALHSPG HVIQPSQALA AGFHRGKQED VHEFLMFTVD AMKKAACLP GH KQVDHHCKDT TLIHQIFGGC WRSQIKCLHC HGISDTFDPY LDIALDIQAA QSVKQALEQL VKPEELNGEN AYHCGLCLQR APASNTLTLH TSAKVLILVL KRFSDVAGNK LAKNVQYPEC LDMQPYMSQQ NTGPLYVYVLY AVL VHAGWSC HDGHYFSYVK AQEVQWYKMD DAEVTVCSSII SVLSQQAYVL FYIQKSEWER HSESVSRGRE PRALGAEDTD RRAKQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTPEF NVGKVEGTL PNALVIHQSK YKCGMKNHHP EQQSSLLNLS STTRTDQESM NTGTLASLQG RTRRAKGKNK HSKRALLVCQ
<b>Source</b>	E.coli
<b>Target Names</b>	USP17L1
<b>Protein Names</b>	Recommended name: Putative ubiquitin carboxyl-terminal hydrolase 17-like protein 1 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 17-like protein 1 Ubiquitin thioesterase 17-like protein 1 Ubiquitin-specific-processing p
<b>Expression Region</b>	1-530
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.