



Recombinant Human Putative ubiquitin carboxyl-terminal hydrolase 17-like protein 1 (USP17L1P)

Product Code	CSB-MP755915HU
Abbreviation	USP17L1P
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7RTZ2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGDDSLYLGG EWQFNHFSKL TSSRPDAAFA EIQRTSLPEK SPLSSETRVD LCDDLAPVAR QLAPREKLPL SSRRPAAVGA GLQNMGNTCY ENASLQCLTY TLPLANYMLS REHSQTCQRP KCCMLCTMQA HITWALHSPG HVIQPSQALA AGFHRGKQED VHEFLMFTVD AMKKAACLP GH KQVDHHCKDT TLIHQIFGGC WRSQIKCLHC HGISDTFDPY LDIALDIQAA QSVKQALEQL VKPEELNGEN AYHCGLCLQR APASNTLTLH TSAKVLILVL KRFSDVAGNK LAKNVQYPEC LDMQPYMSQQ NTGPLYVYVLY AVLHAGWSC HDGHYFSYVK AQEVQWYKMD DAEVTVCSII SVLSQQAYVL FYIQKSEWER HSESVSRGRE PRALGAEDTD RRAKQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTPEF NVGKVEGTL PNALVIHQSK YKCGMKNHHP EQQSSLLNLS STTRTDQESM NTGTLASLQG RTRRAKGKNK HSKRALLVCQ
Source	Mammalian cell
Target Names	USP17L1
Protein Names	Recommended name: Putative ubiquitin carboxyl-terminal hydrolase 17-like protein 1 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 17-like protein 1 Ubiquitin thioesterase 17-like protein 1 Ubiquitin-specific-processing p
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.