



Recombinant Danio rerio UPF0609 protein C4orf27 homolog (zgc:101819)

Product Code	CSB-BP755969DIL
Abbreviation	zgc:101819
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SXS8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAGRGKRKPR SLPQTETPNG EVKKAKEGLK DDKTSVGEEM REEVERLYKL RMPDDFFQFW DFCEGLNADC PQDALKNTLG LQLVGPFDIL SKKHKNSSSQ PNFHLHWRYF YDPPEFQTII QGNADTQHMM GYFRDLPDAL PVFIGENEAK KGYTITQLGD NIFAAVLLFL QKKKKEKRQQ KDDAALNRLE EDLKREAERL GLPLEQKTKS MKQRERKVVT KTFHGAGIVV PVDKNDVGYR ELPESDASLK KICKAIAEAK DDEERMKAFA PIQEMITFVQ FANDECDYGM GYELGIDLFC YGSHYFFKVV RQLLPMAYNL LKRGLEGEIL EAHLASRSQD NLDQLAAV
Source	Baculovirus
Target Names	hpf1
Protein Names	Recommended name: UPF0609 protein C4orf27 homolog
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.