



Recombinant Danio rerio Transcription factor Sox-21-A (sox21a)

Product Code	CSB-MP755981DIL
Abbreviation	sox21a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SZS1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAKPMDHVKR PMNAFMVWSR AQRRKMALDN PKMHNSEISK RLGGEWKLLS DSEKRPFIDE AKRLRAVHMK EHPDYKYRPR RKPKNLIKGD RYHFNVTYNL GEGDPLKSAR LSGDALSDSL SAEKTSVAAA ATRMFFSHHL SANPYPFDL NSKISELPPA PFPHYSVLGY PSGLPAPFGM GVLGAPHAH PSAGSPGHML PYNCLGWPGS GPAVPSSYVL LPGMTKAQHE PYLAYAATI
Source	Mammalian cell
Target Names	sox21a
Protein Names	Recommended name: Transcription factor Sox-21-A Alternative name(s): SRY-box containing gene 21a
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.