



Recombinant Danio rerio Cysteine and histidine-rich domain-containing protein 1 (chordc1)

Product Code	CSB-YP755996DIL
Abbreviation	chordc1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7T3F7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSLMCYNKGC GQRFDLEKNS DDACTYHPGV PVFHDALKGW SCCKRRTTDF SDFLSIAGCT KGPHNQEKPA EPVKPEVKSS LDKNDVKPKY DEFIIQAPKP LESIQRSPD EPFSILQSKI SPSLKQALEK LKLTEENAHE IKEEDSDEIK IGTSCKNGGC CKTFAGQASD DETCLYHAGV PIFHEGMKYW SCCKRKTSDF NTFLSQEGCT RGKHQWKKKD AGKKVVPCRF DWHQTGSQVI ISYAKNSVP ELSLVEGNST VLKIHIFEG EKEFEKQISL WGVIDPSKSL VNMMATKIEI VLKKAEPMSW ARLDLPPVAP PKEKEKEKDV DSEDECDDDE D
Source	Yeast
Target Names	chordc1
Protein Names	Recommended name: Cysteine and histidine-rich domain-containing protein 1 Alternative name(s): CHORD domain-containing protein 1 Morgana
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.