



Recombinant Rat Splicing factor U2AF 26 kDa subunit (U2af114)

Product Code	CSB-EP756019RA-B
Abbreviation	U2af114
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TP17
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	AEYLASIFG TEKDKVNCSF YFKIGACRHG DRCSRLHNKP TFSQTIVLLN LYRNPQNTAQ TADGSHCHVS DVEVQEHYDN FFEEVFTTELQ EKYGEIEEMN VCDNLGDHLV GNVYVKFRRE EDAERAVAEL NNRWFNGQAV HAELSPVTFD RESCCRQYEM GECTRGGFCN FMHLRPISRN LRRQLYGRGP RHRSPPRSHT GHRPRERNRR RSPDHRHGRF
Source	E.coli
Target Names	U2af114
Protein Names	Recommended name: Splicing factor U2AF 26 kDa subunit Alternative name(s): Liver regeneration-related protein LRRG157/LRRG158 U2 auxiliary factor 26 U2 small nuclear RNA auxiliary factor 1-like protein 4
Expression Region	2-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.