



Recombinant Mouse Exosome complex component RRP40 (Exosc3)

Product Code	CSB-BP756031MO
Abbreviation	Exosc3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TQK4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AEVLSAGPE SVAGCRARAV HKVLNQVVLP GEELVLPDHE DVDGLGGAGE QPLRLNAGAR PRLRVVCGPG LRRCGDRLLV TKCGRLRHKE PSGGGGGVYW VDSQQKRYVP VKGDHVGIV IAKSGDIFKV DVGGSEPASL SYLAFEGATK RNRPNVQVGD LIYGQCVVAN KDMPEMVC I DSCGRANGMG VIGQDGLL FK VTLGLIRKLL APDCEIVQEL GKLYPLEIVF GMNGRIWVKA KTIQQTLILA NVLEACEHMT TEQRKQIFAR LAES
Source	Baculovirus
Target Names	Exosc3
Protein Names	Recommended name: Exosome complex component RRP40 EC= 3.1.13.- Alternative name(s): Exosome component 3 Ribosomal RNA-processing protein 40
Expression Region	2-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.