



Recombinant Mouse Interferon kappa (Ifnk)

Product Code	CSB-EP756038MO-B
Abbreviation	Ifnk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TSL0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	LNCVYLDDS ILENVKLLGS TMTGFPLRCL KDITDFKFPK EILPYIQHMK REINAVSYRI SSLALTIFNL KGSIPPVTEE HWERIRSGLF KQVRQAQECF MDEEKENREH PHSEDFLTVY LELGKYFFRI KKFLINKKYS FCAWKIVTVE IRRCFIIFSK SRKLLKMISE SPTFKQELK
Source	E.coli
Target Names	Ifnk
Protein Names	Recommended name: Interferon kappa Short name= IFN-kappa
Expression Region	22-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the type I interferon family. Type I interferons are a group of related glycoproteins that play an important role in host defenses against viral infections. This protein is expressed in keratinocytes and the gene is found on chromosome 9, adjacent to the type I interferon cluster.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.