



# Recombinant L-arabinose isomerase (araA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP756139SZB   |
| <b>Abbreviation</b>      | araA  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q7UDT4  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Shigella flexneri   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MAIFDNYEVW FVIGSQHLYG PETLRQVTQH AEHVVNALNT EAKLPYKLVL<br>KPLGTTPEI TAICRDANYD DRCAGLVVWL HTFSPAKMWI NGLTMLNKPL<br>LQFHTQFNAALPWDSIDMDF MNLNQTAHGG REFGFIGARM RQQHAVVTGH<br>WQDKQAHERI GSWMRQAVSK QDTRHLKVCR FGDNMREVAV<br>TDGDKVAAQI KFGFSVNTWA VGYLVQVVNS ISDGDVNALV DEYESCYTMT<br>PATQIHGEKR QNVLEAARIE LGMKRFLEQG GFHAFTTTFE DLHGLKQLPG<br>LPVQRLMQQG YGFAGEGDWK TAALLRIMKV MSTGLQGGTS FMEDYTYHFE<br>KGNDLVLGSH MLEVCPPIAV EEKPIDVQV LGIGGKDDPA RLIFNTQTGP<br>AIVASLIDLG DRYRLLVNCI DTVKTPHSLP KLPVANALWK AQPDLPTASE<br>AWILAGGAHH TVFSHALNLN DMRQFAEMHD IEITVIDNDT RLPFAKDALAR<br>WNEVYYGFRR |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | araA  |
| <b>Protein Names</b>     | Recommended name: L-arabinose isomerase EC= 5.3.1.4   |
| <b>Expression Region</b> | 1-500   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |