



Recombinant Rat Neurogenin-1 (Neurog1)

Product Code	CSB-MP015730RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70595
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MPAPLETCLS DLDCASSNSG SDLSSFLTDE EDCARLQPLA STSGLSVPAR RSAPTLSGAS NVPGGQDEEQ ERRRRRGRAR VRSEALLHSL RRSRRVKAND RERNRMHNLN AALDALRSVL PSFPDDTKLT KIETLRFAYN YIWALAETLR LADQGLPGGG ARERLLPPQC VPCLPGPPSP ASDTESWGSG AAASPCATVA SPLSDPSSPS ASEDFTYGGP GPLFSFPGLP KDLLHTTPCF IPYH
Source	Mammalian cell
Target Names	Neurog1
Protein Names	Recommended name: Neurogenin-1 Short name= NGN-1 Alternative name(s): Neurogenic basic-helix-loop-helix protein Neurogenic differentiation factor 3 Short name= NeuroD3
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.