



# Recombinant Danio rerio Mediator of RNA polymerase II transcription subunit 6 (med6)

<b>Product Code</b>	CSB-EP757109DIL
<b>Abbreviation</b>	med6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6IQ63
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATVDLRDNL LGISWVDSGW VPILNPSNVL EYFSERSNPF YDRTCNNEVV KMQRSTLDHL TQLVGVEYIL LHTQEPILYI IRKQQRQSPT QVIPLADYYI IAGVVYQAPD LGSVISSRAL SAVHGIQSAF DEAMSFERYH PSKGYWWHFK DQEEKERVKP KSKRKEEPSS LFQRQRVDTL LLDLRNKFPP TFYQTKPGEK PVPVEVKKEP EVPVETVKPQ EERETKPPAP PAPP RPPPTQ TPNKPPPEKR ARVQ
<b>Source</b>	E.coli
<b>Target Names</b>	med6
<b>Protein Names</b>	Recommended name: Mediator of RNA polymerase II transcription subunit 6 Alternative name(s): Mediator complex subunit 6
<b>Expression Region</b>	1-254
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.