



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor ERF086 (ERF086)

Product Code	CSB-MP757125DOA
Abbreviation	ERF086
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6J9Q2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSTSKTLDHN KPFETSQTQM GFALIHQNTS ANTTTTTTTGG ERRGRRSKQA EPGRFLGVRR RPWGRYAAEI RDPTTKERHW LGTFDTAHEA ALAYDRAALS MRGTQARTNF VYTPTDVHTI LTNPNLHSLI VSPYNNNQSF LPNSSPQFVI DHHPHYQNYH QPQQPKHTLP QTVLPAASFK TPVRHQSVDI QAFGNSPQNS SSNGSLSSSL DEENNFFFSL TSEEHNKSNNSG YLDCIVP NHCLKPPPEA TTTQNGAGAS FTTPVASKAS EPYGGFSNSY FEDGEMMMMN HHEFGSCDLS AMITNYGAAA ASMSMEDYGM MEPQDLSSSS IAAFQDVVAD TTGFYSVF
Source	Mammalian cell
Target Names	ERF086
Protein Names	Recommended name: Ethylene-responsive transcription factor ERF086
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.