



# Recombinant Human Protein FAM53A (FAM53A)

<b>Product Code</b>	CSB-EP757499HU-B
<b>Abbreviation</b>	FAM53A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NSI3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVTLITEKLQ SQSLDDLTK AEAGPLQYSA ETLNKSGRLF PLELNDQSPW KVFSGGPPVR SQAATGPDFS FLPGLSAAAH TMGLQWQPQS PRPGAGLGAA STVDPSESTG SSTAPPTKRH CRSLSEPEEL VRCRSPWRPG SSKVWTPVSK RRCDSSGGSAT RQGSPGAVLP RSAVWSTGPT SPATPRPSSA SGGFVDSSEG SAGSGPLWCS AESCLPSTRR RPSLSQERLA GAGTPLPWAS SSPTSTPALG GRRGLLRCS QPCVLSGKRS RRKRRREEDA RWTRPSLDLFL KMTQTLKNSK SLCSLNYEDD DEDDTPVKTV LSSPCDSRGL PGITMPGCSQ RGLRTSPVHP NLWASRESVT SDGSRRSSGD PRDGDSVGEE GVFPRARWEL DLEQIENN
<b>Source</b>	E.coli
<b>Target Names</b>	FAM53A
<b>Protein Names</b>	Recommended name: Protein FAM53A Alternative name(s): Dorsal neural-tube nuclear protein
<b>Expression Region</b>	1-398
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.