



Recombinant *Xenopus laevis* Nucleotide exchange factor SIL1 (sil1)

Product Code	CSB-MP757511XBE
Abbreviation	sil1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NUA7
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	DRSPEYALT KVTEAEDDGN VEAIVENEPD PEDQEIFYPT REWQVVKPGQ AVPAGLHMRL NLQTGKNEAK LLEEKEGKDR PKQKRKSSTT SDKYTKKELK EALTKFKEGA EQLSPAEEKD YLQDIKQFRF PIEDLQKAFN DLNINVETDF EIMTKIVNRF NSSSSTTEKV SALYDLEYVYV HQVDNAQNLL KLGALQLLIN SLNSTDTLLI ENSAFVIGSA LSSNPKVQIE AFEAGALQKL LVILAADQEV SVKKKTLYAL SSMLRQFPYA QQRFMKLGGL QILKNFFKEK NAESLYIRVI TLLYDMIMEK MLLYKENNTE QYEQKYQQYN QINLLESITE QGWCPIISDL LRLPENSRE KVLKAVLTLI PLCRAEFLKD CNLLTLLNSL RKEYEGLAAE ELRSGEQDGY FSEILSLTIS ITDKLK
Source	Mammalian cell
Target Names	sil1
Protein Names	Recommended name: Nucleotide exchange factor SIL1
Expression Region	32-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.