



Recombinant Danio rerio Ribosome biogenesis protein wdr12 (wdr12)

Product Code	CSB-EP757525DIL-B
Abbreviation	wdr12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NX08
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSQLQARFFT ENKKYLVDDV PFSIPASSEV ADLSGVINKL LEAKNAGHKH VEFDLVRGQ FLRTSLDKHI ETENISTEEV VEIEYVEKFT APQPEECIMH DDWISSVEAD SEWILTGSYD KTARIWSMEG KAVMTVAGHA DVVKDVAWVK RDGLNSVLLT ASLDQTVMLW EWNSEFNKVK ARHCCRGHAG SVDTIAVDPT RTKFCSGSWD KMLKIWSAVP TEEDEIEEP NRPRKKQKTE QMGLTRTPLM TLSGHNEAVS SVLWMDADEL CSASWDHTIR LWDAETGGQK STLSGSKVFN CISYSPLCRR LASGSTDRHV RLWDPRSKDG SLVLLSLTSH NGWVTAVRWA PSHEHQLVSG SLDNVVKLWD TRSCKAPLYD LAAHEDKVLK ADWTENGLIL SGGADNKLYT YRYAACVSDA GA
Source	E.coli
Target Names	wdr12
Protein Names	Recommended name: Ribosome biogenesis protein wdr12 Alternative name(s): WD repeat-containing protein 12
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.