



Recombinant Danio rerio F-box/SPRY domain-containing protein 1 (fbxo45)

Product Code	CSB-EP757537DIL
Abbreviation	fbxo45
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NZ03
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSGAVATGGA QSAGLGAAAA GCCSSASGAG SAALLSGGSG IAGRLPSRVL EHIFSYLELM DLLNCSLVCW HWNNCLSDEN SEVWRSLCAR SLSDEAMRSD ILCNLASYKG KLKSFQHALS SHDCSRNVYI KKNGFTHLRN PIAQSTDGAR GKIGFSEGRH AWEIWWEGPL GTVAVIGIAT KRAPMQCQGY VALLGSDDQS WGWNLVDNNL LHNGEVNGNF PQCNNAPKYQ IGERIRVILD MDDKTLA FER GFEFLGVAFR GLPKTCLFPA VSAVYGNTTEV TMVYLGRPLD G
Source	E.coli
Target Names	fbxo45
Protein Names	Recommended name: F-box/SPRY domain-containing protein 1 Alternative name(s): F-box only protein 45
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.