



Recombinant Human F-box and leucine-rich protein 22 (FBXL22)

Product Code	CSB-EP757543HU-B
Abbreviation	FBXL22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P050
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWPLLLTMHIT QLNRECLLHL FSFLDKDSRK SLARTCSQLH DVFEDPALWS LLHFRSLTEL QKDNFLLGPA LRSLICWHS SRVQVCSIED WLKSAFQRSI CSRHESLVND FLLRVCDRLS AVRSPRRREA PAPSSGTPIA VGPKSPRWGG PDHSEFADLR SGVTGARAAA RRGLGSLRAE RPSETPPAPG VSWGPPPPGA PVVISVKQEE GKQGR TGRRS HRAAPPCGFA RTRVCPPTFP GADAFPQ
Source	E.coli
Target Names	FBXL22
Protein Names	Recommended name: F-box/LRR-repeat protein 22 Alternative name(s): F-box and leucine-rich repeat protein 22
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.