



Recombinant Danio rerio Nucleoporin-like protein 2 (nupl2)

Product Code	CSB-MP757548DIL
Abbreviation	nupl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P0U9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MPVCNFFLQG RCRYGDTCWN EHPTGGRGGD YRGNQQPSNR GGFGNRVWIN PSQRGGGGSS SAGGSNEWGR GAASARDVQS SEFSFSQNR SALETQRAGA EDTHHTLDTI QKEMEVWQTS GQWPFSCYSA VNRQISGFIE LCPEELRLEY YTSRASGDIQ PYINSVQQLA NQWRSRVQEL RNMSSSTQIS VIAELKSSSP PASAPGFGSP GPGFGSATSG FGNTSLSSPP AGFGGAGFGS GPQSSSTFSF AQSKTDFGAS NTQQASGFGS SAFAQPSSGF GNPAPSAASF SFAAADSESK PSAGGFGSAS GFSFKTATAG QGSGFGSGFG SSSGFGSSSG FGSSSGFGSA FGSAAPAQSS SFGSTGGAAD TQSGHGLFTA NSELTPEELK EFMAKRFTLG QIPLRPPPAD LLMI
Source	Mammalian cell
Target Names	nupl2
Protein Names	Recommended name: Nucleoporin-like protein 2
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.