



Recombinant Human General transcription factor IIH subunit 2-like protein (GTF2H2C)

Product Code	CSB-MP757553HU
Abbreviation	GTF2H2C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P1K8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDEEPTKR WEGGYERTWE ILKEDESGSL KATIEDILFK AKRKRVEFHH GQVRLGMMRH LYVVVDGSRT MEDQDLKPNR LTCTLKLEY FVEEYFDQNP ISQIGIIVTK SKRAEKLTEL SGNPRKHITS LKEAVDMTCH GEPSLYNSLS MAMQTLKHMP GHTSREVLII FSSLTTC DPS NIYDLIKTLK AAKIRVSVIG LSAEVRVCTV LARETGGTYH VILDESHYKE LLTHHLSPPP ASSSSECSLI RMGFPQHTIA SLSDQDAKPS FSMALHDGNT EPGLTLGGYF CPQCRAKYCE LPVECKICGL TLVSAPHLAR SYHHLFPLDA FQEIPLEEYN GERFCYGCQG ELKDQHVYVC AVCQNVFCVD CDVHVHDSLH CCPGCIHKIP APSGV
Source	Mammalian cell
Target Names	GTF2H2C;
Protein Names	Recommended name: General transcription factor IIH subunit 2-like protein Alternative name(s): General transcription factor IIH polypeptide 2-like protein
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.