



# Recombinant Human Centrosomal protein of 85 kDa (CEP85)

<b>Product Code</b>	CSB-EP757559HU
<b>Abbreviation</b>	CEP85
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6P2H3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAMQEKYPTE GISHVTSPSS DVIQKGSSLG TEWQTPVISE PFRSRF SRCS SVADSGDTAI GTSCSDIAED FCSSSGSPPF QPIKSHVTIP TAHVMPSTLG TSPA KPNSTP VGPSSSKLPL SGLAESVGMT RINGDLGAMKH SPGLSRDLMY FSGATGENGI EQSWFPAVGH ERQEEARKFD IPSMESTLNQ SAMMETLYSD PHHRVRFHNP RTSTSKEL YR VLPEAKKAPG SGAVFERNGP HSNSSGVLPL GLQPAPGLSK PLPSQVWQPS PDTWHPREQS CELSTCRQQL ELIRLQMEQM QLQNGAICHH PAAFGPSLPI LEPAQWISIL NSNEHLLKEK ELLIDKQRKH ISQLEQKVRE SELQVHSALL GRPAPFGDVC LLRLQELQRE NTFLRAQFAQ KTEALSREKI DLEKKLSASE VEVQLIRESL KVALQKHSEE VKKQEERVKG RDKHINNLKK KCQKESEQNR EKQQRIETLE RYLADLPTLE DHQKQSQQLK DSELKSTELQ EKVTELESLL EETQAICREK EIQLES LRQR EAEFSSAGHS LQDKQSVEET SGEGPEVEME SWQKRYDSLQ KIVEKQQQKM DQLRSQVQSL EQEVAQE EGT SQALREEAQR RDSALQQLRT AVKELSVQNQ D LIEKNLTLQ EHLRQAQPGS PPSPTAQLA LELHQELASC LQDLQAVCSI VTQRAQGHDP NLSLLLGIHS AQHPETQLDL QKPDVIKRL EEVQQLRRDI EDLRTTMSDR YAQDMGENCV TQ
<b>Source</b>	E.coli
<b>Target Names</b>	CEP85
<b>Protein Names</b>	Recommended name: Centrosomal protein of 85 kDa Alternative name(s): Coiled-coil domain-containing protein 21
<b>Expression Region</b>	1-762
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.