



Recombinant Human DENN domain-containing protein 1B (DENND1B)

Product Code	CSB-EP757568HU
Abbreviation	DENND1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P3S1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDCRTKANPD RTFDLVLKVK CHASENEDPV VLWKFPEDFG DQEILQSVPK FCFPFDVERV SQNQVGQHFT FVLTDIESKQ RFGFCRLTSG GTICLCILSY LPWFEVYYKL LNTLADYLAK ELENLNETL RSLYNHPVVK ANTPVNLSVN QEIFIACEQV LKDQPALVPH SYFIAPDVTG LPTIPESRNL TEYFVAVDVN NMLQLYASML HERRIVISS KLSTLTACIH GSAALLYPMY WQHIYIPVLP PHLLDYCCAP MPYLIGIHSS LIERVKNKSL EDVVMLNVDN NTLESPFSDL NNLPSDVVSA LKNKLLKQST ATGDGVARAF LRAQAALFGS YRDALRYKPG EPITFCEESF VKHRSSVMKQ FLETAINLQL FKQFIDGRLA KLNAGRGFSD VFEEEITSGG FCGGKDKLQY DYPFSQ
Source	E.coli
Target Names	DENND1B
Protein Names	Recommended name: DENN domain-containing protein 1B Alternative name(s): Connecden 2 Protein FAM31B
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.