



Recombinant Mouse SH3KBP1-binding protein 1 (Shkbp1)

| | |
|--------------------------|--|
| Product Code | CSB-EP757592MO |
| Abbreviation | Shkbp1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6P7W2 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AVATTAVEG VPSRGAPGEV IHLNVGGKRF STSRQTLTWI PDSFFSSLLS GRISTLKDET GAIFIDRDPT VFAPILNFLR TKELDPRGVH GSSLLHEAQF YGLTPLVRRL QVREELDRSS CGNVLFNGYL PPPVFPVKRR NRHSLVGPQQ IGGRPAPVRR SNTMPPNLGN AGLLGRMLDE RAPPSPSGQP EEPGMVRLVC GHHNWIAVAY THFLVCYRLK EASGWQLAFS SPRLDWPIER LALAARVLGG APGEHDKMVA AATGSEILLW ALQAQGGGSE IGVFHLGVPV EALFFVGNQL IATSHTGRIG VWNAVTKHWQ VQEVQPITSY DAAGSFLLLG CSNGSIYYVD VQKFPLRMKD NDLLVSELYR DPAEDGVTAL SVYLTPTKSD SGNWIEIAYG TSSGVVRVIV QHPETVGS GP QLFQTF SVHR SPVTKIMLSE KHLISVCADN NHVRTWSVTR FRGMISTQPG STPLASF KIL ALESVDGLGG CSAGNDIGPY GERDDQQVFI QKVVPNASQL FVRLSSTGQR VCSVRSVDGS PTTAFTVLEC EGSRRLGSRP RRYLLTGQAN GSLAMWDLTT AMDGLGQTPA GGLTEELMD QLEQCELSPL ASSRGSFPSP SPRTSLTSLH SASSNTSLYG PRGSPSPQA EARRRGAGSF VDRCQELARG APELRWPPTP APRPSTSLGN PLTALKKTLN ETSF |
| Source | E.coli |
| Target Names | Shkbp1 |
| Protein Names | Recommended name: SH3KBP1-binding protein 1 Alternative name(s): SETA-binding protein 1 |
| Expression Region | 2-704 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.