



Recombinant Mouse SH3KBP1-binding protein 1 (Shkbp1)

Product Code	CSB-MP757592MO
Abbreviation	Shkbp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P7W2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AVATTAVEG VPSRGAPGEV IHLNVGGKRF STSRQTLTWI PDSFFSSLLS GRISTLKDET GAIFIDRDPT VFAPILNFLR TKELDPRGVH GSSLLHEAQF YGLTPLVRRL QVREELDRSS CGNVLFNGYL PPPVFPVKRR NRHSLVGPQQ IGGRPAPVRR SNTMPPNLGN AGLLGRMLDE RAPPSPSGQP EEPGMVRLVC GHHNWIAVAY THFLVCYRLK EASGWQLAFS SPRLDWPIER LALAARVLGG APGEHDKMVA AATGSEILLW ALQAQGGGSE IGVFHLGVPV EALFFVGNQL IATSHTGRIG VWNAVTKHWQ VQEVQPITSY DAAGSFLLLG CSNGSIYYVD VQKFPLRMKD NDLLVSELYR DPAEDGVTAL SVYLTPTKSD SGNWIEIAYG TSSGVVRVIV QHPETVGS GP QLFQTF SVHR SPVTKIMLSE KHLISVCADN NHVRTWSVTR FRGMISTQPG STPLASFKIL ALESVDGLGG CSAGNDIGPY GERDDQQVFI QKVVPNASQL FVRLSSTGQR VCSVRSVDGS PTTAFTVLEC EGSRRLGSRP RRYLLTGQAN GSLAMWDLTT AMDGLGQTPA GGLTEELMD QLEQCELSPL ASSRGSFPSP SPRTSLTSLH SASSNTSLYG PRGSPSPQA EARRRGAGSF VDRCQELARG APELRWPPTP APRPSTSLGN PLTALKKTLN ETSF
Source	Mammalian cell
Target Names	Shkbp1
Protein Names	Recommended name: SH3KBP1-binding protein 1 Alternative name(s): SETA-binding protein 1
Expression Region	2-704
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.