



# Recombinant Human Connector enhancer of kinase suppressor of ras 3 (CNKSR3)

<b>Product Code</b>	CSB-EP757602HU
<b>Abbreviation</b>	CNKSR3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6P9H4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPVTKWSPK QVVDWTRGLD DCLQQYVHKF EREKINGEQL LQISHQDLEE LGVTRIGHQE LVLEAVDLLC ALNYGLETDN MKNLVLKLRA SSHNLQNYIS SRRKSPAYDG NTSRKAPNEF LTSVVELIGA AKALLAWLDR APFTGITDFS VTKNKIIQLC LDLTTTVQKD CFVAEMEDKV LTVVKVLNGI CDKTIRSTTD PVMSQCACLE EVHLPNIKPG EGLGMYIKST YDGLHVITGT TENSPADRSQ KIHAGDEVIQ VNQQTVVGWQ LKNLVKKLRE NPTGVVLLK KRPTGSFNFT PAPLKNLRWK PPLVQTSPPP ATTQSPESTM DTSLKKEKSA ILDLYIPPPP AVPYSPRDEN GSFVYGGSSK CKQPLPGPKG SESPNSFLDQ ESRRRRFTIA DSDQLPGYSV ETNILPTKMR EKTPSYGKPR PLSMPADGNW MGIVDPFARP RGHGRKGEDA LCRYFSNERI PPIIESSSP PYRFSRPTTE RHLVRGADYI RGSRCYINS D LHSSATIPFQ EEGTKKKS GS SATKSSSTEP SLLVSWFTRL KLLTH
<b>Source</b>	E.coli
<b>Target Names</b>	CNKSR3
<b>Protein Names</b>	Recommended name: Connector enhancer of kinase suppressor of ras 3 Short name= Connector enhancer of KSR 3 Alternative name(s): CNK homolog protein 3 Short name= CNK3 CNKSR family member 3 Maguin-like protein
<b>Expression Region</b>	1-555
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.