



Recombinant Danio rerio Src kinase-associated phosphoprotein 2 (skap2)

Product Code	CSB-MP757644DIL
Abbreviation	skap2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PG29
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MRSVPEDFTT LISDLETFLT DVLKGENLSK KAKEKKDAFI KRIKDVKTSY AQEFKDRRDE EFPEPDEADT NDGGLHSEQ TDRDDENAYD GVQQSPPVAA QDLQAVLKSG YLEKRRKDHS FFGNEWQKRW CALNNSIFY YGSEKDKQQK GEFNIVGYTV KMNNTLRKDA KRDCCEVSA PDKRVYQFCA ASEKEAKEWV EHIDFLIKDL GGIPEDEEE YDDCLSMSQS EVAILPDDDI YEELPEEDVP QPSKPKVTLV NKPPPPATPV TAAVNKSTDY ANYFQGLWDC IGDQPDELSF KRGDIYILS KEYDMFGWWV GEMKGVIGIV PKEYLLELYV L
Source	Mammalian cell
Target Names	skap2
Protein Names	Recommended name: Src kinase-associated phosphoprotein 2 Alternative name(s): Src family-associated phosphoprotein 2
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.