



Recombinant Human BTB/POZ domain-containing protein KCTD18 (KCTD18)

Product Code	CSB-BP757654HU
Abbreviation	KCTD18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PI47
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGHKAEDEV LDVLRNLVGG CIYTARRESL CRFKDSMLAS MFSGRFPLKT DESGACVIDR DGRLFKYLLD YLHGEVQIPT DEQTRIALQE EADYFGIPYP YLSLDHLANE METYSLRSNI ELKKALDFC DSYGLVCNKP TVWVLHYLNT SGASCESRII GVIYATKTDGT DAIEKQLGGR IHSKGIFKRE AGNNVQYIWS YYSVAELKKM MDAFDAWEGK GVSYWRVPHE LIECWTLER PLLGSLRHMA PIRKRRITF NEADESVNYK TGPKPVRFLG PSTSTQIKVK NSASVTVSPA SAIQTSAGAT ANRFQSGSRR KAAQRSAPSR ATALVGTGAP GHPQASPGAA SAENGGTHLP PAKVLLSDKK PTPQRVIKLN RTPLCATAPC LPSPTATRQA NSLKPLPGEA ARALGVRTEN GKNKGN
Source	Baculovirus
Target Names	KCTD18
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD18
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.