



Recombinant Human Forkhead box protein R2 (FOXR2)

Product Code	CSB-BP757660HU
Abbreviation	FOXR2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PJQ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDLKLKDCEF WYSLHGQVPG LLDWDMRNEL FLPCTTDQCS LAEQILAKYR VGVMKPPEMP QKRRPSPDGD GPPCEPNLWM WVDPNILCPL GSQEAPKPSG KEDLTNISPF PQPPQKDEGS NCSEDKVVES LPSSSSEQSP LQKQGIHSPS DFELTEEEAE EPDDNSLQSP EMKCYQSQKL WQINNQEKS W QRPLNCSHL IALALRNNPH CGLSVQEIYN FTRQHFFFFW TAPDGWKSTI HYNLCFLDSF EKVPDSLKDE DNARPRSCLW KLTKEGHRRF WEETRVLAFA QRERIQECMS QPELLTSLFD L
Source	Baculovirus
Target Names	FOXR2
Protein Names	Recommended name: Forkhead box protein R2 Alternative name(s): Forkhead box protein N6
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.