



Recombinant Human E3 ubiquitin-protein ligase TRAF7 (TRAF7)

Product Code	CSB-YP757671HU
Abbreviation	TRAF7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Q0C0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSGKSARYN RFSGGPSNLP TPDVTTGTRM ETTFGPAFSA VTTITKADGT STYKQHCRT P SSSSTLAYSP RDEEDSM PPI STPRRSDSAI SVRSLHSESS MSLRSTFSLP EEEEEPEPLV FAEQPSVKLC CQLCCSVFKD PVITTCGHTF CRRCALKSEK CPVDNVKLT V VNNIAVAEQ IGELFIHCRH GCRVAGSGKP PIFEVDPRGC PFTIKLSARK DHEGSCDYRP VRCPPNNPSCP PLLRMNLEAH LKECEHIKCP HSKYGCTFIG NQDTYETHLE TCRFEGLEKEF LQQTDDRFHE MHVALAQKDQ EIAFLRSM LG KLSEKIDQLE KSLELKFDVL DENQSKLSED LMEFRRDASM LNDELSHINA RLNMGILGSY DPQQIFKCKG TFGVGHQGPVW CLCVYSMGDL LFSGSSDKTI KVWDTCTTYK CQKTLEGH DG IVLALCIQGC KLYSGSADCT IIVWDIQNLQ KVNTIRAH DN PVCTLVSSH N VLFSGSLKAI KVWDIVGTEL KKKKELTGLN HWVRALVAAQ SYLYSGSYQT IKIWDIRTL D CIHVLQTSGG SVYSIAVTNH HIVCGTYENL IHVWDIESKE QVRTLTGHVG TVYALAVIST PDQTKVFSAS YDRSLRVWSM DNMICTQTLL RHQGSVTALA VSRGRLFSGA VDSTVKVWTC
Source	Yeast
Target Names	TRAF7
Protein Names	Recommended name: E3 ubiquitin-protein ligase TRAF7 EC= 6.3.2.- Alternative name(s): RING finger and WD repeat-containing protein 1 RING finger protein 119 TNF receptor-associated factor 7
Expression Region	1-670
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.