



Recombinant *Saccharomyces cerevisiae* Transposon Ty4-H Gag polyprotein (TY4A-H)

Product Code	CSB-EP757678SVG-B
Abbreviation	TY4A-H
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Q5P6
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MATPVRDETR NVIDDNISAR IQSKVKTNDT VRQTPSSLRK VSIKDEQVKQ YQRNLNRFKT ILNGLKAEEL KLSETDDIQM LAEKLLKLGE TIDKVENRIV DLVEKIQLLE TNENNNILHE HIDATGTYYL FDTLTSTNKR FYPKDCVFDY RTNNVENIPI LLNFKKFIK KYQFDDVFEN DIIIDPREN EILCKIIKEG LGESLDIMNT NTTDIFRIID GLKNKYRSLH GRDVRIRAWK KVLVDTTCRN SALLMNKLQK LVLMEKWIFS KCCQDCPNLK DYLQEAIMGT LHESLRNSVK QRLYNIPHNV GINHEEFLIN TVIETVIDLS PIADDQIENS CMYCKSVFHC SINCKKKPNR ELRPDSTNFS KTYYLQGAQR QQQLKSSAKR TKVLEQDTTK VKQSVQQQKT GNY
Source	E.coli
Target Names	TY4A-H
Protein Names	Recommended name: Transposon Ty4-H Gag polyprotein Short name= TY4A Short name= Transposon Ty4 protein A
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.