



Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 6 (USP17L6)

Product Code	CSB-YP757690HU
Abbreviation	USP17L6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6QN14
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDDSL YLRG EWQFNHFSKL TSSRPDAAFA EIQR TSLPEK SPLSCETRVD LCDDLAPVAR QLAPREKLPL SSRPAAVGA GLQNMGNTCY VNASLQCLTY TPPLAN YMLS REHSQTCHRH KGCMLCTMQA HITRALHNPG HVIQPSQALA AGFHRGKQED AHEFLMFTVD AMKKA CLPGH KQVDHHSKDT TLIHQIFGGY WRSQIKCLHC HGISDTFDPY LDIALDIQAA QSVQQALEQL VKPEELNGEN AYHCGVCLQR APASKTLTLH TSAKVLILVL KRFS DVTGNK IAKNVQYPEC LDMQP YMSQQ NTGPLVYVLY AVL VHAGWSC HNGHYFSYVK AQEGQWYKMD DAEVTASSIT SVLSQQAYVL FYIQKSEWER HSESVSRGRE PRALGSED
Source	Yeast
Target Names	USP17L6P
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase 17-like protein 6 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 17-like protein 6 Ubiquitin thioesterase 17-like protein 6 Ubiquitin-specific-processing protease 1
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.