



Recombinant Human POTE ankyrin domain family member G (POTEG)

Product Code	CSB-EP757729HU-B
Abbreviation	POTEG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6S5H5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVAEAGSMPA ASSVKKPFGL RSKMGKWCRH CFPWCRGSGK SNVGTSGDHD DSAMKTLRSK MGKWCRHCFP WCRGSSKSNV GTSGDHDDSA MKTLRSKMGK WCCHCFPCCR GSGKSKVGPW GDYDDSAFME PRYHVRREDL DKLHRAAWWG KVPRKDLIVM LKDTDMNKKD KQKRTALHLA SANGNSEVVK LLLDRRCQLN ILDNKKRTAL TKAVQCREDE CALMLLEHGT DPNIPDEYGN TALHYAIYNE DKLMKALLL YGADIESKNK HGLTPLLLGV HEQKQQVVKF LIKKKANLNA LDRYGRTALI LAVCCGSASI VSLLEQNID VSSQDLGQT AREYAVSSH NVICQLLSYD KEKQMLKVSS ENSNPEQDLK LTSEEEESQRL KGSENSQPEE MSQEPEINKG GDRKVEEEMK KHGSTHMGFP ENLPNGATAD NGDDGLIPPR KSRTPESSQQF PDTENEQYHS DEQNNTQKQL SEEQNTGILQ DEILIHEEKQ IEVAENEF
Source	E.coli
Target Names	POTEG
Protein Names	Recommended name: POTE ankyrin domain family member G Alternative name(s): ANKRD26-like family C member 2 Prostate, ovary, testis-expressed protein on chromosome 14 Short name= POTE-14
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.