



Recombinant Rat NF-kappa-B essential modulator (Ikbkg)

| | |
|--------------------------|--|
| Product Code | CSB-MP757759RA |
| Abbreviation | Ikbkg |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6TMG5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Rattus norvegicus (Rat) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSRHLWKNQL SEMVQPSGGP AEDQDMLGEE SSLGKPAMLH LPSEQGTPET LQRCLEENQE LRDAIRQSNQ MLRERCEELL HFQVSQREEK EFLMCKFQEA RKLVERLSLE KLDLRRQREQ ALEDLEHLKK CQQQMAEDKA SVKAQVTSLL GELQESQSRL EAATKERQTL EGRIRAVSEQ VRQLESEREV LQQQHSVQVD QLRMQNQSVE AALRMERQAA SEEKRKLAQL QAAYHQLFQD YDSHIKSSKG MQLEDLRQQL QQAEEALVAK QELIDKLKEE AEQHKIVMET VPVLKAQADI YKADFQAERH AREKLVERKE LLQEQLEQLQ REFNKLKVGVC HESARIEDMR KRHVETSQPP LLPAPAHHSF HLALSNQRRS PPEEPPDFCC PKCQYQAPDM DTLQIHVMEC IE |
| Source | Mammalian cell |
| Target Names | Ikbkg |
| Protein Names | Recommended name: NF-kappa-B essential modulator Short name= NEMO Alternative name(s): Ikb kinase-associated protein 1 Short name= IKKAP1 Inhibitor of nuclear factor kappa-B kinase subunit gamma Short name= I-kappa-B kinas |
| Expression Region | 1-412 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.