



Recombinant Human CMRF35-like molecule 5 (CD300LD), partial

Product Code	CSB-MP757797HU1
Abbreviation	CD300LD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6UXZ3
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AK ITGPTTVNGS EQGSLTVQCA YGSGWETYLK WRCQGADWNY CNILVKTNGS EQEVKKNRVS IRDNQKNHVF TVTMENLKRD DADSYWCGTE RPGIDLGVKV QVTINPGTQT AVSEWTTTTA SLAFTAAATQ KTSSPLTRSP LKSTH
Source	Mammalian cell
Target Names	CD300LD
Protein Names	Recommended name: CMRF35-like molecule 4 Short name= CLM-4Alternative name(s): CD300 antigen-like family member D CMRF35-A4 CD_antigen= CD300d
Expression Region	19-165
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.