



# Recombinant *Naumovozyma castellii* Dihydroorotate dehydrogenase (fumarate) (URA1)

<b>Product Code</b>	CSB-BP757802NAY
<b>Abbreviation</b>	URA1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6V3X0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Naumovozyma castellii</i> (strain ATCC 76901 / CBS 4309 / NBRC 1992 / NRRL Y-12630) (Yeast) ( <i>Saccharomyces castellii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPASLTTKFL NQTYENPLMN ASGVHCMTTK ELDELAASRA GAFITKSSTT LERAGNPEPR YVSVPLGSIN SMGLPNQGID YYLEYVLERQ RTHPNEPAIF FSVAGMSIDE NINMLRKIQD SEFHGITELN LSCPNVPGKP QVAYDFDLTK EILDKVFQFF KKPLGIKLP YDFAHFDIM AGILNQYPLA YVNSINSVGN GLYIDVDKES VVKPKNGFG GIGGEYVKPT ALANVRAFYT RLNPSIKIVG TGGIRTGKDV FEHLLCGASM VQIGTELVKE GVPIFERLER ELKEVMDKKG YTTIEEFRGK LNSL
<b>Source</b>	Baculovirus
<b>Target Names</b>	URA1
<b>Protein Names</b>	Recommended name: Dihydroorotate dehydrogenase (fumarate) Short name= DHOD Short name= DHODase Short name= DHODEHASE EC= 1.3.98.1 Alternative name(s): Dihydroorotate oxidase
<b>Expression Region</b>	1-314
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.