



Recombinant Human Transcription factor LBX2 (LBX2)

Product Code	CSB-MP757859HU
Abbreviation	LBX2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6XYB7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNSGREPRTP RTLLSIADIL APRMVPRAPS APQLPESGPG PTSPLCALEE LTSKTFRGLD ARALQPSEGR AGPDALGPGP FGRKRRKSRT AFTAQQVLEL ERRFVFQKYL APSERDGLAT RLGLANAQVV TWFQNRRAKL KRDVEEMRAD VASLRALSPE VLCSLALPEG APDPGLCLGP AGPDSRPHLS DEEIQVDD
Source	Mammalian cell
Target Names	LBX2
Protein Names	Recommended name: Transcription factor LBX2 Alternative name(s): Ladybird homeobox protein homolog 2
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.