



# Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 15 (Os02g0567200, LOC\_Os02g35910)

<b>Product Code</b>	CSB-EP757895OFG
<b>Abbreviation</b>	Os02g0567200, LOC_Os02g35910
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6YTI2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTRSKSVPV PAGGGAATVP LAVLLRREVV SEKTAAERPE LQVGLFSQAK KGEDYTFLKP DCERLPGVPS SSFSAFGLFD GHNGNGAAIY TKENLLSNIL TAIPADLNRE DWLAALPRAM VAAFVKTDKD FQTKARSSGT TVTFVIIDGL FITVASVGDS RCVLEAEGSI YHLSADHRFD ASKEEVDRVT ESGGDVGRLN VVGGAIEIGPL RCWPGGLCLS RSIGDQDVGQ FIVVPYVKQ VKLSTAGGRL IISSDGVWDV LTAEVAFNCS RTLPEAAAE QIVKEAVQQK GLRDDTTCIV VDILPDKANL TMPHTKKQPG MGVFKNMFRK KTPSDSSSHT DREYMDPDIV EEIFEDGCAF LSKRLDSEYP VRNMFKLFIC AICQVELKPS QGISVHEDSS QPGNLRRWDG PFLCQGCQEK KEAMEGKRRS RDSSSRNSGS SE
<b>Source</b>	E.coli
<b>Target Names</b>	Os02g0567200
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 15 Short name= OsPP2C15 EC= 3.1.3.16
<b>Expression Region</b>	1-442
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.