



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 15 (Os02g0567200, LOC_Os02g35910)

Product Code	CSB-MP757895OFG
Abbreviation	Os02g0567200, LOC_Os02g35910
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YTI2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSTRSKSVPV PAGGGAATVP LAVLLRREVV SEKTAAERPE LQVGLFSQAK KGEDYTFLKP DCERLPGVPS SSFSAFGLFD GHNGNGAAIY TKENLLSNIL TAIPADLNRE DWLAALPRAM VAAFVKTDKD FQTKARSSGT TVTFVIIDGL FITVASVGDS RCVLEAEGSI YHLSADHRFD ASKEEVDRVT ESGGDVGRLN VVGGAIEIGPL RCWPGGLCLS RSIGDQDVGQ FIVVPYVKQ VKLSTAGGRL IISSDGWVDV LTAEVAFNCS RTLPEAAAE QIVKEAVQQK GLRDDTTCIV VDILPDKANL TMPHTKKQPG MGVFKNMFRK KTPSDSSSHT DREYMDPDIV EEIFEDGCAF LSKRLDSEYP VRNMFKLFIC AICQVELKPS QGISVHEDSS QPGNLRRWDG PFLCQGCQEK KEAMEGKRRS RDSSSRNSGS SE
Source	Mammalian cell
Target Names	Os02g0567200
Protein Names	Recommended name: Probable protein phosphatase 2C 15 Short name= OsPP2C15 EC= 3.1.3.16
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.