



Recombinant *Oryza sativa* subsp. japonica Clathrin light chain 3 (Os07g0461500, LOC_Os07g27750)

Product Code	CSB-EP757914OFG-B
Abbreviation	Os07g0461500, LOC_Os07g27750
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z3A8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSSFDSVA AV AGDGDADDDD VLPPAPFDPA ADGAQGGLGA LRRGHRFATS YSSFGTAASE DDLAGAGAGT DGGVGAGIPL GSSSNGGAAY GYGGSGDVMN GHVDQIGDVM GGGVVVGDGG GIDDDL FAGA GDGDDGPVLP PPEAMKEEGI LRREWRRQNA LMLEEKERKE RERRGEIIAE ADEFKRSFAE KRKLNGDTNR AQNRDREKLF LAKQEKFHGE AEKQYWKAIA EMVPHEIPGL EKRGRREKQ SAEANAKAKQ PGVVVVQGTK PGKPTDL SRM RQVLMKLKQT PPPHMAPP PP QPAKDTGGDT DANKDGEAEK AAGEIEK KAA GGEKEAAAGP PVTAAAAADA QANKAAAEET AKK
Source	E.coli
Target Names	Os07g0461500
Protein Names	Recommended name: Clathrin light chain 3
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.