



# Recombinant *Oryza sativa* subsp. japonica Endoglucanase 4 (GLU14)

<b>Product Code</b>	CSB-EP757919OFG
<b>Abbreviation</b>	GLU14
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6Z715
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FNYADALD KAILFFEAQR SGKLPQGQRV AWRADSGLSLSD GSADGVDLGAG GYYDAGDNVK FGLPMAFTVT MLSWSVIEFG DMMPARRSSF LGGIFGGGGV AQLDNARAASV RWGADYLLKA ATATPDTLYV QVADPYQDHR CWERPEDMDT PRSVYKVTPQ SPGSDVAGET AAALAAASIV FRVSDPSYSA KLLDAAQLVF DFADKYRGSY SDSLSSVVCV FYCSHSYHDE LLWAASWLHL ASPEKKDVYL SYIGSNGHAL GAEQDDFTFS WDDKRVATKG FLQSRADGLQ LYKAHTDNYI CSLVPGANGF QSQYTPGGLL FKEGDSNMQY VTSTAFLLLT YAKYLSSSAA TVSCGSTAVS PSTLISLAKK QVDYILGANP AGMSYVMVGF ARYPRHVHHR GASMPVSRDH PARIGCDEGF RYLHSPEPDR NLLAGAVVGG PDAGDAFADG RDNYAQAEPS TYTNAPLVGA LAFFAGAHKI FTP
<b>Source</b>	E.coli
<b>Target Names</b>	GLU14
<b>Protein Names</b>	Recommended name: Endoglucanase 4 EC= 3.2.1.4 Alternative name(s): Endo-1,4-beta glucanase 4 OsCel9F OsGLU14
<b>Expression Region</b>	23-503
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.