



Recombinant *Oryza sativa* subsp. japonica Probable uridine nucleosidase 1 (URH1)

Product Code	CSB-YP757928OFG
Abbreviation	URH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZJ05
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGSNEQIHRD KLIIDTPGI DDSMTILMAF RAPTVEIIGL TTIFGNTTTK NATQNALLLC ERAGHPEVPV AEGSAEPLKG GEPRVADFVH GSDGLGNLFL PAPTSKKVDE NAAEFMVNKV SQFPGEVSIL ALGPLTNVAL AIKRDPFAS KVKKIVVLGG AFFAAGNVSP AAANIYGDP EAADIVFTSG ADVDVVGINI TTQVCFTDED LLELRNSK GK HAQFLCDMCQ FYRDWHAESD GFHGIFLHDP VSFTALVHPE YFTFKKGVVR VETQGICTGH TLMDQGLKKW NSENPWSGYK PISVAWTVDV PNVLAFVKEL LMAP
Source	Yeast
Target Names	URH1
Protein Names	Recommended name: Probable uridine nucleosidase 1 EC= 3.2.2.3 Alternative name(s): Uridine ribohydrolase 1
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.